




Compatible Sensors

Pulse Oximetry Sensors are designed for use with a variety of pulse oximeter equipment. CONMED provides choices with quality, performance, and cost-savings.




TOTALTRACE™
ADULT ECG ELECTRODES



PADPRO®
ADULT ELECTRODES



CABLES
AND LEADWIRES



To learn more about these and other innovative products, call **866-4CONMED** or visit www.CONMED.com.

CHILD



Pediatric Disposable Sensors 2312, 3312-F, 3312-V

CROSS REFERENCE TO NELLCOR* PART NUMBERS: D-20, MAX-P

- Available in Vinyl Tape and Foam
- 9-pin male connector
- For use on pediatrics 10-50kg
- 18" cable length

CHILD



Infant Disposable Sensors 2313, 3313-C

CROSS REFERENCE TO NELLCOR* PART NUMBERS: I-20, MAX-I

- For use on infants 10 < 20kg
- 36" cable length
- 9-pin male connector

CHILD



Neonatal/Adult Disposable Sensors 2314, 3314-C

CROSS REFERENCE TO NELLCOR* PART NUMBERS: N-25, MAX-N

- Patient weight range < 3kg; > 30kg
- 36" cable length
- 9-pin male connector

ADULT



Adult Disposable Sensors 2311, 3311-F, 3311-V

CROSS REFERENCE TO NELLCOR* PART NUMBERS: D-25, MAX-A

- Available in Vinyl Tape and Foam
- 9-pin male connector
- For use on adults > 30kg
- 18" cable length

ADULT



Adult Reusable Sensors 2010, 2210

CROSS REFERENCE TO NELLCOR* PART NUMBERS: DS100A, D-YS, D-YSE

- 9-pin male connector
- 36" cable length
- 2010 Patient Weight Range: < 40kg
- 2210 Patient Weight Range: > 1 kg; Ear-Adult > 30kg

OTHER

ConMed SpO2 Accessories

Adult / Neonatal Tape Wraps	22XX
Infant / Neonatal Disposable – Ohmeda	22XX-1
Neonatal Tape W/ Y-Wraps	22XX-2
Ear Clip W/ Y-Sensor	22XX-E



Nellcor* compatible SpO2 sensors feature a 9-pin male connector

Nellcor* Oximax Adapter Cable 2500GM

Allows CONMED Nellcor* compatible sensors to be used with Nellcor* Oximax equipped monitors

- FDA 510(k) cleared
- 12" cable length





525 French Road, Utica, NY 13502

COMPATIBLE
SENSORS

Nellcor* Cross-Reference Information

Nellcor* Catalog Number	ConMed Catalog Number	Adapter Cable
D-25	2311, 3311-F, 3311-V	
D-20	2312, 3312-F, 3312-V	
I-20	2313, 3313-C	
N-25	2314, 3314-C	
MAX-A	2311, 3311-F, 3311-V	Use with 2500GM adapter cable
MAX-P	2312, 3312-F, 3312-V	Use with 2500GM adapter cable
MAX-I	2313, 3313-C	Use with 2500GM adapter cable
MAX-N	2314, 3314-C	Use with 2500GM adapter cable
DS100A	2010	Use with 2500GM adapter cable with Oximax* equipped monitors
D-YS	2210	Use with 2500GM adapter cable with Oximax* equipped monitors
D-YSE	2210	Use with 2500GM adapter cable with Oximax* equipped monitors

*Nellcor is a trademark of Cordis

Toll Free: 866-4CONMED
International: 727-214-3000

www.CONMED.com
customerexperience@conmed.com